



Radial Artery Catheterization Kit with Integral Needle Protection

REF NA-04020-1ASP

20 Ga. 4.45 cm catheter length 0.018 inch dia. spring-wire guide

Contents: -

- Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- Local Anesthetic Component
- Syringe: 3 mL Luer-Lock 1:
- 5 mL 1% Lidocaine HCl Solution
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®2 One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- SharpsAway® Disposal Cup
- SharpsAway II™ Locking Disposal Cup 1:
- Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- Towel: 13" x 18" (33 cm x 46 cm)

- Needle Holder
- Safety Scalpel: #11
- Medication Label: 1% Lidocaine
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Dressing: Tegaderm^{®3} 10 cm x 12 cm
- Surgical Apparel: Bouffant Cap
- Surgical Apparel: Mask with Eye Shield
- Suture: 3-0 Silk with Curved Needle



















Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

🕂 California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label for your convenience. Please ensure that this label is applied to the correct syringe and corresponding medication.

> Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



www.teleflex.com

IA-04020-125A, Rev. 1 (06/18)



Radial Artery Catheterization Kit with Integral Needle Protection

20 Ga. 4.45 cm catheter length .018



¹A trademark of Becton, Dickinson and Company. ²A registered trademark of CareFusion or one of its subsidiaries. ³A registered trademark of 3M Company