



ARROW

Radial Arterial Catheterization Kit with Sharps Safety Features

REF NA-04020-K

20 Ga. | 4.45 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
- 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Wing Clip: Suture
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Towel
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex




Store between 15 - 30°C (59 - 86°F). Avoid freezing and excessive heat above 40°C (104°F)

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.



Manufacturer: Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Teleflex®

www.teleflex.com

IA-04020-139C, Rev. 01 (2020-12)



ARROW

Radial Arterial Catheterization Kit with Sharps Safety Features

REF NA-04020-K

20 Ga. | 4.45 cm catheter length

Contents:

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
- 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: Wing Clip: Suture
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Towel
- 1: Drape: 18" x 26" with 2" x 3" fenestration, with adhesive
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)

<sup>1</sup> A trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.



Caution



Consult instructions for use



Do not reuse



Do not resterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry



Do not use if package is damaged



Not made with natural rubber latex




Store between 15 - 30°C (59 - 86°F). Avoid freezing and excessive heat above 40°C (104°F)

Rx only. Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Fluid path components are non-pyrogenic.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.



Manufacturer: Arrow International LLC Subsidiary of Teleflex Incorporated 3015 Carrington Mill Blvd. Morrisville, NC 27560 USA

Teleflex®

www.teleflex.com

IA-04020-139C, Rev. 01 (2020-12)



Radial Arterial Catheterization Kit with Sharps Safety Features

REF NA-04020-K

20 Ga. | 4.45 cm catheter length

ARROW



Radial Arterial Catheterization Kit with Sharps Safety Features

REF NA-04020-K

20 Ga. | 4.45 cm catheter length

ARROW