



# Arterial Catheterization Kit

REF NA-04020-1A

20<sup>Ga.</sup> | 4.45 cm catheter length | .018 inch dia. spring-wire guide

**Contents:**

- 1: Catheter: 20 Ga. x 4.45 cm Radiopaque Polyurethane over 22 Ga. TW Introducer Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: Wing Clip: Suture
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Needle Holder

- 1: Safety Scalpel: #11
- 1: Medication Label: 1% Lidocaine
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: Tegaderm®<sup>3</sup> 10 cm x 12 cm
- 1: Surgical Apparel: Bouffant Cap
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Suture: 3-0 Silk with Curved Needle

<sup>1</sup> A trademark of Becton, Dickinson and Company.  
<sup>2</sup> A registered trademark of CareFusion or one of its subsidiaries.  
<sup>3</sup> A registered trademark of 3M Company



STERILE EO



~~LATEX~~  
 Not made with natural rubber latex

20°C (68°F) - 25°C (77°F)  
 Store between 20 - 25°C (68 - 77°F).

Rx only.

**Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**California Prop. 65**

**Warning: Cancer and Reproductive Harm [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).**

Fluid path components are non-pyrogenic.

Arrow International provides the enclosed medication label for your convenience. Please ensure that this label is applied to the correct syringe and corresponding medication.

Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IA-04020-140H (06/18)



# Arterial Catheterization Kit

REF NA-04020-1A

20<sup>Ga.</sup> | 4.45 cm catheter length | .018 inch dia. spring-wire guide