

REF NA-04220-S1A

ARROW



# QuickFlash® Radial Artery Catheterization Kit

20<sup>Ga.</sup> | 3.81 cm catheter length | .018 inch dia. spring-wire guide

## Contents:

- 1: QuickFlash® Radial Artery Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane over 21 Ga. TW Needle with Integral .018" (0.46 mm) dia. Spring-Wire Guide
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: SharpsAway® Disposal Cup
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 36" x 41" (91 cm x 104 cm) with 3" x 5" (8 cm x 13 cm) fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Checklist Sheet
- 1: Medication Label: 1% Lidocaine
- 1: Dressing: GRIP-LOK®<sup>3</sup> Catheter Fixation Device
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: I.V. Advanced 10 cm x 12 cm
- 1: Dressing: BIOPATCH®<sup>4</sup>
- 1: Surgical Apparel: Mask with Eye Shield
- 1: Surgical Apparel: Bouffant Cap
- 1: Tubing<sup>2</sup>: 6" Extension with One-Way Stopcock

<sup>3</sup>A registered trademark of TIDI.

<sup>4</sup>A registered trademark of Johnson & Johnson Corporation.

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



**LATEX**  
Not made with natural rubber latex

**PHT DEHP**  
Precaution: Contains Phthalate: DEHP

20°C (68°F) | 25°C (77°F)  
Store between 20-25°C (68-77°F).

REF NA-04220-S1A

ARROW



# QuickFlash® Radial Artery Catheterization Kit

20<sup>Ga.</sup> | 3.81 cm catheter length | .018 inch dia. spring-wire guide

Manufacturer:  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IA-04220-125A, Rev. 01 (2020-12)