

REF NA-04550-K

ARROW



Arterial Catheterization Kit

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque Polyurethane with Integral Suture Wing
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) RW
- 1: Filter Needle: 18 Ga. TW x 1-1/2" (3.81 cm) with Blunt Fill Tip
- 1: Syringe: 3 mL Luer-Lock
- 2: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®² One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Towel: 13" x 18" (33 cm x 46 cm)
- 1: Scalpel: #11
- 4: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)

¹A trademark of Becton, Dickinson and Company.

²A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



STERILE EO



~~LATEX~~
Not made with natural rubber latex

Store between 15°C (59°F) and 30°C (86°F). Avoid freezing and excessive heat above 40°C (104°F).

REF NA-04550-K

ARROW



Arterial Catheterization Kit

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Manufacturer:
Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex
www.teleflex.com

IA-04550-105A, Rev. 01 (2020-12)