

REF CA-09803

# Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port  
for use with 7 - 7.5 Fr. Catheters

# ARROW®

**Contents:**

- One: Sheath: 8.5 Fr. x 4" (10 cm)  
Arrow-Flex® Radiopaque Polyurethane  
with Integral Hemostasis Valve/Side Port,  
and Tissue Dilator
  - One: Spring-Wire Guide, Marked: .035"  
(.89 mm) dia. x 17-13/16" (45 cm)  
(Straight Soft Tip on One End – "J" Tip on Other)  
with Arrow Advancer™<sup>1</sup>
  - One: Catheter: 18 Ga. x 2-1/2" (6.35 cm)  
Radiopaque over 20 Ga. RW Introducer Needle
  - One: Pressure Transduction Probe
  - One: Introducer Needle: 18 Ga. x 2-1/2"  
(6.35 cm) XTW and 5 mL Arrow® Raulerson  
Spring-Wire Introduction Syringe
  - One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)  
and 5 mL Luer-Slip Syringe
  - One: Injection Needle: 25 Ga. x 1" (2.54 cm)  
and 3 mL Luer-Lock Syringe
  - One: Suture: 3-0 Braided Silk with Straight Needle
  - One: Dust Cap: Non-Vented
  - One: Cath-Gard® Catheter Contamination Shield:  
80 cm with TwistLock™ Adapter
  - One: Drape: 24" x 36" with 4" fenestration
  - One: Drape: 24" x 36"
  - One: CSR Wrap
  - Two: Gauze Pads: 2" x 2"
  - Five: Gauze Pads: 4" x 4"
  - One: Obturator: 8 Fr.
  - One: Scalpel: #11
  - One: Stopcock: 4-way High-Flow
  - One: SharpsAway® Disposal Cup
- <sup>1</sup> US Patent Nos. 5,484,419.

**Rx only.**


**Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains dry natural rubber.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA



IA-09803-101D (8/09)



REF CA-09803    8.5 Fr. 10 cm sheath length    .035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit

Peel to Open