

# Two-Lumen CVC Kit

with Blue FlexTip®, ARROWgãrd Blue® Catheter



## Arrow Select Kits

# ASK-22702-CA

> ARROWgãrd Blue® Catheter

### Contents:

- 1: Two-Lumen Indwelling Catheter: 7 Fr. x 6-5/16" (16 cm) Radiopaque Polyurethane with Blue FlexTip®, ARROWgãrd® Antimicrobial Surface Treatment<sup>†</sup>, Extension Line Clamps, Injection Site Caps**
- 1:** Spring-Wire Guide, Marked: .025" (0.64 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1:** Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1:** Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) RW
- 1:** Injection Needle: Eclipse™<sup>2</sup> 25 Ga. x 1" (2.54 cm) RW
- 1:** Pressure Transduction Probe
- 1:** Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1:** Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW and 5 mL Luer-Slip Syringe
- 1:** Tissue Dilator
- 1:** SecondSite™ Adjustable Hub: Catheter Clamp and Fastener
- 1:** SharpsAway II™ Locking Disposal Cup
- 1:** CSR Wrap
- 1:** Needle Holder
- 1:** Safety Scalpel: #11
- 2:** Gauze Pad: 2" x 2"
- 5:** Gauze Pad: 4" x 4"
- 1:** Suture: 3-0 Silk with Curved Needle

<sup>1</sup> Licensed under US Patent No. 6,706,024.

<sup>2</sup> A trademark of Becton, Dickinson and Company.

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**Contraindications: The ARROWgãrd Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**


Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic.

**Store below 25°C ( 77°F). Avoid excessive heat above 40°C ( 104°F).**



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA  
www.arrowintl.com

IA-22702-102A (5/15)



Catheter Cross Section

Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (18 Ga.)	0.33	1300
Proximal (14 Ga.)	0.50	5600

\*Priming volumes are approximate and are done without an injection site cap. Injection site cap priming volume is 0.17 mL.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

## Two-Lumen CVC Kit

with Blue FlexTip®, ARROWgãrd Blue® Catheter

ASK-22702-CA

**7** Fr. **2** Lumen **16** cm catheter length **.025** inch diameter spring-wire guide

> ARROWgãrd Blue® Catheter

Arrow Select Kits

