



Central Venous Catheterization Kit with Blue FlexTip® Arrowg+ard Blue® Catheter

REF CA-24301

16^{Ga.} 20^{cm} catheter length

Contents:

- Indwelling Catheter: 16 Ga. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment¹, Integral Extension Line, Extension Line Clamp, Integral Suture Wing
- Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- **Pressure Transduction Probe** 1:
- Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe 1:
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

- Tissue Dilator: 7 Fr. (2.3 mm) x 10.2 cm
- Dust Cap: Non-Vented 1:
- SecondSite[™] Adjustable Hub: Fastener 1:
- SecondSite™ Adjustable Hub: Catheter Clamp 1:
- SharpsAway® Disposal Cup
- Drape: 21" \times 36" (53 cm \times 91 cm) with 3" (7.6 cm) fenestration, with adhesive
- Scalpel: #11 1:
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent No. 6,706,024.




















Not made with natural rubber latex



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).



Lumen	Priming Volume* (mL)	Flow Rate [†] (mL/hr)	
Distal (16 Ga.)	0.43	2781	
	0.15		_

Priming volumes are approximate and are done without accessories.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Fluid path components are non-pyrogenic.

Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

> Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IA-24301-201B (01/20)

ARROW





Central Venous Catheterization Kit with Blue FlexTip® Arrowg+ard Blue® Catheter

catheter

Flow rate values are approximate and are determined using deionized water at 100 cm head