



ARROWg⁺ard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

REF **CA-25122-F**

12 Fr. | **2** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- 1: Two-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWg⁺ard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 2: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent No. 6,706,024.



Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**


California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: The Arrow® large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.

The ARROWg⁺ard Blue® antimicrobial surface catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



| Lumen | Priming Volume* (mL) | Flow Rate† (mL/min) | Venous Pressure Setting (mmHg) |
|-------------------|----------------------|---------------------|--------------------------------|
| Distal (12 Ga.) | 1.60 | 200 | 100 |
| | | 340 | 200 |
| Proximal (12 Ga.) | 1.50 | 200 | 100 |
| | | 340 | 200 |

* Priming volumes are approximate and are done without accessories.
 † Blood flow rates are approximate and were determined by varying the venous pressure setting from 100 to 200 mmHg in an in vitro experiment using blood analog solution.

DuPont™ Tyvek® | made in USA

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IA-25122-101D (6/17)



ARROWg⁺ard Blue® Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

REF **CA-25122-F**

12 Fr. | **2** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

