



Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter for High Volume Infusions

REF **CA-25123-F** **12**^{Fr.} | **3**^{Lumen} | **20**^{cm} catheter length | **.035**^{inch dia.} spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, ARROWgard® Antimicrobial Surface Treatment¹, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Tissue Dilator: Stepped
- 3: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹Licensed under US Patent No. 6,706,024.

DuPont™ Tyvek® packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Sterile. Fluid path components are non-pyrogenic.

Contraindications: The ARROWgard Blue® antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine acetate, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).



Caution



Consult instructions for use



Do not reuse



Do not sterilize



Sterilized by ethylene oxide



Keep away from sunlight



Keep dry




Do not use if package is damaged



Not made with natural rubber latex



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/min.)
Distal (16 Ga.)	0.40	49
Medial (12 Ga.)	1.30	196
Proximal (12 Ga.)	1.20	221

* Priming volumes are approximate and are done without accessories.

† Based on 100 cm head height.

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IA-25123-101E (4/17)



Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® ARROWgard Blue® Catheter for High Volume Infusions

REF **CA-25123-F** **12**^{Fr.} | **3**^{Lumen} | **20**^{cm} catheter length | **.035**^{inch dia.} spring-wire guide