Product No. SB-09883-SPercutaneous Sheath
Introducer Kit
with Integral Hemostasis Valve/Side Port
for use with 7.5 Fr. CathetersARRR/COLSPANE/Side Port
for use with 7.5 Fr. Catheters

Special Order Product

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex[®] Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 1-3/4" (4.45 cm) Radiopaque over 20 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 3 mL Luer-Lock Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Luer-Slip Syringe
- 1: Obturator: 8 Fr.
- 1: Cath-Gard[®] Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
 - Dust Cap: Non-Vented
 - SharpsAway® Disposal Cup
 - Drape: 24" x 36" with 4" fenestration
- 1: Drape: 24 1: CSR Wrap

1:

1:

1:

1:

5:

1:

- Drape: 21" x 36" with 3" fenestration, with adhesive
- 1: Scalpel: #11
 - Stopcock: 4-way High-Flow
- 2: Gauze Pad: 2" x 2"
 - Gauze Pad: 4" x 4"
 - Suture: 3-0 Braided Silk with Straight Needle

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use - do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.

This Product Contains Dry Natural Rubber.





Arrow International, Inc. 2400 Bernville Road Reading, PA 19605 USA www.teleflex.com

