

PSI Kit

with Integral Hemostasis Valve/Side Port
(Recommended for use with 7 Fr. Catheters)



ARROW
INTERNATIONAL

Arrow Select Kits

ASK-09907-SB

> Sharps Safety Features

> Latex-Free

Contents:

- 1: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .025" (.64 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- 1: Pressure Transduction Probe
- 1: **Injection Needle: SafetyGlide®² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe**
- 1: **Injection Needle: SafetyGlide®² 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: Obturator: 8 Fr.
- 1: Dust Cap: Non-Vented
- 1: **SharpsAway IITM³ Locking Disposal Cup**
- 2: CSR Wraps
- 1: Drape: 29" x 18" with adhesive
- 1: **Maximal Barrier Drape™: 54" x 96" with two 4" fenestrations**
- 1: **Safety Scalpel #11**
- 1: Stopcock: 4-way High-Flow
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 2: Prep Sponge Swabs
- 1: Suture: 3-0 Silk with Straight Needle
- 1: Tray: Prep

¹ US Patent No. 5,484,419.

² A registered trademark of Becton, Dickinson and Company.

³ US Patent Nos. 6,123,193 and 6,276,527.

Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway IITM Locking Disposal Cup
- > Safety Scalpel

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

This product contains NO natural rubber latex.

Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA
www.arrowintl.com

IB-09907-101B (8/10)

PSI Kit with Integral Hemostasis Valve/Side Port

ASK-09907-SB

9 Fr.

10 cm
sheath
length

.025 inch diameter
spring-wire
guide

- > Sharps Safety Features
- > Latex-Free

Arrow Select Kits



ARROW
INTERNATIONAL