

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port
for use with 7.5 - 8 Fr. Catheters



ARROW
INTERNATIONAL

Arrow Select Kits

NA-09903-CDC

> **Sharps Safety Features**

Contents:

- 1: Sheath: 9 Fr. x 4-1/8" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Pressure Transduction Probe
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: **Injection Needle: SafetyGlide®¹ 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: **Injection Needle: SafetyGlide®¹ 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe**
- 1: Obturator: 8 Fr.
- 1: 5 mL HCl, 1% Lidocaine Solution
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 1: **3 mL Applicator Pouch, 2% CHG and 70% IPA ChloraPrep®² One-Step Solution with Hi-Lite Orange™ Tint**
- 1: Needle Holder
- 1: **Drape: 41" x 55" with 4" fenestration**
- 1: Tubing: 20" Extension±
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Local Anesthetic Components
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: **Safety Scalpel: #11**
- 1: **SharpsAway II™ Locking Disposal Cup**
- 1: Medication Label: 1% Lidocaine
- 1: Stopcock: 4-way High-Flow

¹ A registered trademark of Becton, Dickinson and Company.

² A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.

Infection Protection Feature:

- > **ChloraPrep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint
- > Extra-Large Drape size

Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway II™ Locking Disposal Cup
- > Safety Scalpel

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

±Precaution: Contains DEHP.

Arrow International provides the enclosed medication label for your convenience. Please ensure that this label is applied to the correct syringe and corresponding medication.

Store at 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).



Arrow International LLC
3015 Carrington Mill Blvd., Morrisville, NC 27560 USA
www.arrowintl.com

IC-09903-104B, Rev. 01 (2020-12)

PSI Kit with Integral Hemostasis Valve/Side Port

NA-09903-CDC

9 Fr.

10 cm sheath length

.035 in. diameter spring-wire guide

> **Sharps Safety Features**

Arrow Select Kits



ARROW
INTERNATIONAL