

MAC™ Two-Lumen Central Venous Access Kit

with Integral Hemostasis Valve for use with 7.5 - 8 Fr. Catheters



ARROW
INTERNATIONAL

Arrow Select Kits

NA-11242-CDC

> **Sharps Safety Features**


Contents:

- 1: Two-Lumen Central Venous Access Device¹: 115 mm Radiopaque Polyurethane with 9 Fr. Distal Lumen, Integral Hemostasis Valve, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Transduction Probe
- 1: Tissue Dilator
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: **Injection Needle: SafetyGlide®² 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe**
- 1: **Injection Needle: SafetyGlide®² 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: **3 mL Applicator Pouch, 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution with Hi-Lite Orange™ Tint**
- 1: 5 mL HCl, 1% Lidocaine Solution
- 1: Needle Holder
- 1: **Drape: 41" x 55" with 4" fenestration**
- 1: Tubing: 20" Extension[±]
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Local Anesthetic Components
- 1: Obturator: 8 Fr.
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Medication Label: 1% Lidocaine
- 1: **SharpsAway II™ Locking Disposal Cup**
- 1: SharpsAway® Disposal Cup
- 1: **Safety Scalpel: #11**
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- 2: Dust Caps: Non-Vented
- 1: Dust Cap: Vented
- 1: Stopcock: 4-way High-Flow

¹ Licensed under US Patent No. 5,250,038.

² A registered trademark of Becton, Dickinson and Company.

³ A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)
Distal (9 Fr.) no catheter	1.80	30450
Distal (9 Fr.) with 8 Fr. catheter	1.15	12390
Distal (9 Fr.) with 7 Fr. catheter	1.45	19660
Proximal (12 Ga.) no catheter	0.60	11950

* Priming volumes are approximate and are done without an injection cap.
† Based on 100 cm head height.

Infection Protection Features:

- > **ChloraPrep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint
- > Extra-Large Drape size

Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway II™ Locking Disposal Cup
- > Safety Scalpel

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

±Precaution: Contains DEHP.

Arrow International provides the enclosed medication label for your convenience. Please ensure that this label is applied to the correct syringe and corresponding medication.

Store at 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).



Arrow International LLC
3015 Carrington Mill Blvd., Morrisville, NC 27560
www.arrowintl.com

IC-11242-100B, Rev. 01 (2020-12)

MAC™ Two-Lumen Central Venous Access Kit

NA-11242-CDC

9 Fr. distal lumen

2 Lumen

11.5 cm device length

.035 inch diameter spring-wire guide

> **Sharps Safety Features**

Arrow Select Kits



ARROW
INTERNATIONAL