

# Pressure Injectable Multi-Lumen CVC Kit

with Blue FlexTip® Catheter



**ARROW**  
INTERNATIONAL

Arrow Select Kits

## NA-15703-PCDC

> Sharps Safety Features

### Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 8" (20 cm) Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™
- 1: Transduction Probe
- 1: Tissue Dilator
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: **Injection Needle: SafetyGlide®<sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe**
- 1: **Injection Needle: SafetyGlide®<sup>1</sup> 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe**
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: **3 mL Applicator Pouch, 2% CHG and 70% IPA ChloroPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint**
- 1: 5 mL HCl, 1% Lidocaine Solution
- 1: Needle Holder
- 1: **Drape: 41" x 55" with 4" fenestration**
- 1: Tubing: 20" Extension<sup>±</sup>
- 1: Suture: 3-0 Silk with Curved Needle
- 1: Local Anesthetic Components
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 1: Medication Label: 1% Lidocaine
- 1: **SharpsAway II™ Locking Disposal Cup**
- 1: SharpsAway® Disposal Cup
- 1: SecondSite™ Adjustable Hub: Catheter Clamp and Fastener
- 1: **Safety Scalpel: #11**
- 1: Flow Rate Information Card with Injection Log

<sup>1</sup> A registered trademark of Becton, Dickinson and Company.

<sup>2</sup> A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.

### Infection Protection Features:

- > ChloroPrep® One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint
- > Extra-Large Drape size

### Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway II™ Locking Disposal Cup
- > Safety Scalpel

### Rx only.

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**All contents sterile unless otherwise indicated.**

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

**±Precaution: Contains DEHP.**


Arrow International provides the enclosed medication label for your convenience. Please ensure that this label is applied to the correct syringe and corresponding medication.

**Store at 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).**



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA  
www.arrowintl.com

IC-15703-109B (3/16)



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	MAX Pressure Injection Flow Rate‡ (mL/sec)
Distal (16 Ga.)	0.44	3200	10
Medial (18 Ga.)	0.39	1500	5
Proximal (18 Ga.)	0.39	1600	5

\* Priming volumes are approximate and are done without an injection cap.  
† Based on 100 cm head height.  
‡ Pressure injection flow rates are determined using media of 11.8 centipoise viscosity, with 60" administration set tubing.

## Pressure Injectable Multi-Lumen CVC Kit

NA-15703-PCDC

7 Fr.

3 Lumen

20 cm catheter length

.032 inch diameter spring-wire guide

> Sharps Safety Features

Arrow Select Kits



**ARROW**  
INTERNATIONAL