



Radial Artery Catheterization Kit

REF **HF-04020-1-S**

20^{Ga.} | **3.81**^{cm} catheter length | **.025** inch dia. spring-wire guide

Contents:

- | | |
|---|--|
| <ul style="list-style-type: none"> 1: Indwelling Catheter: 20 Ga. x 3.81 cm Radiopaque Polyurethane with Integral Suture Wing 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) 1: Injection Needle: 25 Ga. x 1" (2.54 cm) 1: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW 1: Local Anesthetic Components 1: Syringe: 3 mL Luer-Lock 1: 5 mL 1% Lidocaine HCl Solution | <ul style="list-style-type: none"> 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep^{®1} One-Step Solution 1: SharpsAway II™ Locking Disposal Cup 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration 1: Safety Scalpel: #11 4: Gauze Pad: 2" x 2" (5 cm x 5 cm) 3: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Suture: 3-0 Braided Silk with Straight Needle <p><small>¹A registered trademark of CareFusion or one of its subsidiaries.</small></p> |
|---|--|



Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.
 Fluid path components are non-pyrogenic.

DuPont™ Tyvek®
packaging

Manufacturer:
 Arrow International LLC
 Subsidiary of Teleflex Incorporated
 3015 Carrington Mill Blvd.
 Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IF-04020-105D, Rev. 2 (2020-12)



Radial Artery Catheterization Kit

REF **HF-04020-1-S**

20^{Ga.} | **3.81**^{cm} catheter length | **.025** inch dia. spring-wire guide