

# Arterial Line Kit



Arrow Select Kits

## HF-04510-1

> **Sharps Safety Features**

**Contents:**

- 1: Indwelling Catheter: 20 Ga. x 5" (12 cm) Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm)  
(Straight Soft Tip on One End – "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque  
over 22 Ga. TW Introducer Needle
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>1</sup>  
One-Step Solution
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Suture: 3-0 Braided Silk with Straight Needle
- 1: Local Anesthetic Components
- 3: Gauze Pads: 2" x 2"
- 3: Gauze Pads: 4" x 4"
- 1: **Safety Scalpel: #11**
- 1: **SharpsAway II™ Locking Disposal Cup**

<sup>1</sup> A registered trademark of Carefusion or one of its subsidiaries.



**Rx only.**

**Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**All contents sterile unless otherwise indicated.**

Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic. Contains medication.

**Store between 20 - 25°C (68 - 77°F).**



Arrow International LLC  
3015 Carrington Mill Blvd., Morrisville, NC 27560 USA  
www.arrowintl.com

IF-04510-100E, Rev. 2 (2020-12)

## Arterial Line Kit

HF-04510-1

**20** Ga. **12** cm catheter length **.025** inch diameter spring-wire guide

> **Sharps Safety Features**

Arrow Select Kits

