

REF HF-04510-1

ARROW®



Arterial Line

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Contents:

- 1: Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP
- 1: Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep® One-Step Solution
- 1: SharpsAway® II Locking Disposal Cup
- 1: SharpsAway® Disposal Cup
- 1: Drape: 41" x 55" with 4" fenestration
- 1: Safety Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 3: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.


Fluid path components are non-pyrogenic.



STERILE EO



Not made with natural rubber latex



Store between 20-25°C (68-77°F).

REF HF-04510-1

ARROW®



Arterial Line

20^{Ga.} | 12^{cm} catheter length | .025^{inch dia.} spring-wire guide

Arrow International LLC
Subsidiary of Teleflex Incorporated
3015 Carrington Mill Blvd.
Morrisville, NC 27560 USA

Teleflex®
www.teleflex.com

IF-04510-100F (2021-05)