



ARROWgard Blue[®] Two-Lumen Hemodialysis Catheterization Kit for High Volume Infusions

REF ASK-46122-VG

12^{Fr.} | 2^{Lumen} | 25^{cm} catheter length | .035^{inch dia.} spring-wire guide

Contents:

- | | | |
|---|--|---|
| <ul style="list-style-type: none"> 1: Two-Lumen Indwelling Catheter: 12 Fr. x 10" (25 cm) Radiopaque Polyurethane with Blue FlexTip[®], ARROWgard[®] Antimicrobial Surface Treatment¹, Extension Line Clamps 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1: Pressure Transduction Probe 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe 1: Injection Needle: Eclipse^{™2} 22 Ga. x 1-1/2" (3.81 cm) RW and 5 mL Luer-Slip Syringe 1: Injection Needle: Eclipse^{™2} 25 Ga. x 1" (2.54 cm) RW and 3 mL Luer-Lock Syringe | <ul style="list-style-type: none"> 1: Syringe: 10 mL Luer-Lock 1: Tissue Dilator: 10 Fr. (3.5 mm) x 10.2 cm 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm 1: SharpsAway[®] Disposal Cup 1: SharpsAway II[™] Locking Disposal Cup 1: Maximal Barrier Drape[™] with 4" fenestration 1: Towel 1: Needle Holder 1: Forceps 1: Safety Scalpel: #11 1: Checklist/CLIP Sheet 1: Sterile Procedure Sign 2: Gauze Pad: 2" x 2" (5 cm x 5 cm) | <ul style="list-style-type: none"> 5: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Dressing: Tegaderm^{®3} 10 cm x 12 cm 1: Surgical Apparel: Impervious Gown 1: Surgical Apparel: Bouffant Cap 1: Surgical Apparel: Mask with Eye Shield 1: Suture: 3-0 Braided Silk with Straight Needle 1: Suture: 3-0 Silk with Curved Needle 1: HemoHopper[®] Fluid Receptacle |
|---|--|---|


¹ Licensed under US Patent No. 6,706,024.

² A trademark of Becton, Dickinson and Company.

³ A registered trademark of 3M Company.


Contraindications: The ARROWgard Blue[®] antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).

The Arrow[®] large-bore two-lumen catheter is not designed for long-term hemodialysis or for use in patients with thrombosed vessels.



Lumen	Priming Volume* (mL)
Distal (12 Ga.)	1.4
Proximal (12 Ga.)	1.6

* Priming volumes are approximate and are done without an injection cap.



¹ In vitro tests performed using 3:1 glycerine/saline at 37°C.



Not made with natural rubber latex



Store below 25° C (77° F).
Avoid excessive heat above 40° C (104° F)