

Arterial Line Kit



Arrow Select Kits

ASK-04510-HFH

- > **Sharps Safety Features**
- > **Maximal Barrier Precautions**

Contents:

- 1: Indwelling Catheter: 20 Ga. x 5" (12 cm) Radiopaque FEP
- 1: Spring-Wire Guide, Marked: .025" (.64 mm) x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
- 1: **Injection Needle: SafetyGlide®² 18 Ga. x 1-1/2" (3.81 cm) RW**
- 1: **Injection Needle: SafetyGlide®² 25 Ga. x 1" (2.54 cm)**
- 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW
- 1: Syringe: 3 mL Luer-Lock
- 1: **1 Applicator Pouch, 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution**
- 1: **3 mL Applicator Pouch, 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution with Hi-Lite Orange™ Tint**
- 1: 5 mL HCl, 1% Lidocaine Solution
- 1: **SharpsAway II™⁴ Locking Disposal Cup**
- 1: Local Anesthetic Components
- 1: **Maximal Barrier Drape™: 54" x 96" with 4" fenestration**
- 1: Forceps
- 1: Needle Holder
- 1: **Safety Scalpel: #11**
- 1: Dressing: Tegaderm®⁵ 11.5 cm x 8.9 cm
- 3: Gauze Pads: 2" x 2"
- 3: Gauze Pads: 4" x 4"
- 1: **Surgical Apparel: Bouffant Cap**
- 1: **Surgical Apparel: Impervious Gown**
- 1: **Surgical Apparel: Mask with Eye Shield**
- 1: Suture: 3-0 Silk with Curved Needle

¹ US Patent No. 5,484,419.
² A registered trademark of Becton, Dickinson and Company.
³ A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.
⁴ US Patent Nos. 6,123,193 and 6,276,527.
⁵ A registered trademark of 3M Company.

Infection Protection Features:

- > **ChloroPrep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol
- > **ChloroPrep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint

Maximal Barrier Precautions:

- > Drape Sheet
- > Bouffant Surgical Cap
- > Impervious Surgical Gown
- > Surgical Mask with Eye Shield

Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway II™ Locking Disposal Cup
- > Safety Scalpel

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

Store at 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Arrow International LLC
 3015 Carrington Mill Blvd., Morrisville, NC 27560 USA
 www.arrowintl.com

IH-04510-100D, Rev. 01 (2020-12)

Arterial Line Kit

ASK-04510-HFH

- > **Sharps Safety Features**
- > **Maximal Barrier Precautions**



20 Ga. **12** cm catheter length **.025** inch diameter spring-wire guide