

PSI Kit

with Integral Hemostasis Valve/Side Port
for use with 7 - 7.5 Fr. Catheters



ARROW
INTERNATIONAL

Arrow Select Kits

ASK-09801-UAH

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
- 1: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Pressure Transduction Probe
- 1: Syringe: 3 mL Luer-Lock
- 1: Syringe: 5 mL Luer-Slip
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Cath-Gard® Catheter Contamination Shield: 80 cm with Sterile Tape
- 1: Obturator: 8 Fr.
- 1: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: CSR Wrap
- 1: Drape: 24" x 36"
- 1: Drape: 24" x 36" with 4" fenestration, with adhesive
- 1: Scalpel: #11
- 1: Stopcock: 3-Way High-Flow
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- 2: Prep Sponge Swabs
- 1: Suture: 3-0 Nylon with Straight Needle
- 1: Tray: Prep

¹ US Patent No. 5,484,419.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

This product contains dry natural rubber.

Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA
www.arrowintl.com

IH-09801-102A (5/09)

PSI Kit with Integral Hemostasis Valve/Side Port

ASK-09801-UAH

Arrow Select Kits



ARROW
INTERNATIONAL

8.5 Fr.

10 cm sheath length

.035 inch diameter spring-wire guide