PSI Kit

with Integral Hemostasis Valve/Side Port for use with 7 Fr. Catheters



Arrow Select Kits

ASK-09803-SRH

> Sharps Safety Features

Contents:

- 1: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End − "J" Tip on Other) with Arrow Advancer™
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle 1:
- 1: Pressure Transduction Probe
- 1: Injection Needle: SafetyGlide^{®1} 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Injection Needle: SafetyGlide^{®1} 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow[®] Raulerson Spring-Wire Introduction Syringe
- 1: Cath-Gard[®] Catheter Contamination Shield: 80 cm with TwistLockTM Adapter
- Obturator: 8 Fr. 1:
- Dust Cap: Non-Vented
- 1: SharpsAway IITM Locking Disposal Cup
- CSR Wraps 3:
- Drape: 29" x 18" with adhesive
- 1: Maximal Barrier DrapeTM: 54" x 96" with 4" fenestration
- 1: Safety Scalpel: #11
- Stopcock: 4-way High-Flow
- 2: Gauze Pads: 2" x 2"
- 5: Gauze Pads: 4" x 4"
- Suture: 3-0 Silk Braided with Straight Needle
- $^{\,1}\,$ A registered trademark of Becton, Dickinson and Company.

Safety Features:

- > SafetyGlide® Protected Needles
- > SharpsAway II™ Locking Disposal Cup
- > Safety Scalpel

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication.

This Product Contains Dry Natural Rubber.











Arrow International, Inc. 2400 Bernville Road, Reading, PA 19605 USA www.arrowintl.com

PSI Kit with Integral Hemostasis Valve/Side Port

ASK-09803-SRH

Select

