



# Arrowgard<sup>+</sup> Blue<sup>®</sup> PSI Kit



REF **ASK-29903-HGH**

**9** Fr. | **10** cm | **.035** inch dia.  
sheath length | spring-wire guide

## Contents:

- 1: Sheath: 9 Fr. x 4-1/8" (10 cm) Arrow-Flex<sup>®</sup> Radiopaque Polyurethane with Arrowgard<sup>+</sup> Antimicrobial Surface Treatment<sup>1</sup>, Integral Hemostasis Valve/Side Port and Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Injection Needle: SafetyGlide<sup>™2</sup> 25 Ga. x 1" (2.54 cm)
- 1: Introducer Needle: Echogenic 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow<sup>®</sup> Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: SafetyGlide<sup>™2</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Obturator: 8 Fr.
- 1: Dust Cap: Non-Vented
- 1: SharpsAway<sup>®</sup> Disposal Cup
- 1: SharpsAway II<sup>™</sup> Locking Disposal Cup
- 1: Maximal Barrier Drape<sup>™</sup> with 4" fenestration
- 1: Needle Holder
- 1: Safety Scalpel: #11
- 1: Stopcock: 4-way High-Flow
- 1: Checklist/CLIP Sheet
- 1: Sterile Procedure Sign
- 1: Dressing: Tegaderm<sup>®3</sup> 6 cm x 7 cm
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Silk with Curved Needle
- 1: HemoHopper<sup>®</sup> Fluid Receptacle

<sup>1</sup>Licensed under US Patent No. 6,706,024.

<sup>2</sup>A trademark of Becton, Dickinson and Company.

<sup>3</sup>A registered trademark of 3M Company.

**Contraindications:** The Arrowgard<sup>+</sup> Blue antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).



Not made with natural rubber latex



Store below 25°C (77°F). Avoid excessive heat above 40°C (104°F)