

REF PI-01351-SS

**Latex-Free**

# PICC Seldinger Conversion Set

for use with Arrow® 3 Fr. PICC

# ARROW®

**Contents:**

- One: . . . . . Sheath: 17 Ga. x 1-1/4" (3.18 cm) Radiopaque Peel-Away over 3 Fr. Dilator
- One: . . . . . Spring-Wire Guide: .016" (.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- One: . . . . . Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- One: . . . . . Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW

**Rx only.**


**Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use — do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA



1H01351-101C (4/16)

REF **PI-01351-SS** for use with Arrow **3** Fr. PICC **.016** inch dia. spring-wire guide

**PICC Seldinger Conversion Set**

