

REF PI-01351-LS5

LATEX-FREE

# PICC Seldinger Conversion Set

to be used with Arrow® 3 Fr. PICC

# ARROW®

**Contents:**

- One: .....Sheath: 17 Ga. x 2-3/4" (7 cm) Radiopaque Peel-Away over 3 Fr. Dilator
- One: ..... Spring-Wire Guide: .016" (.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- One: ..... Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- One: .....Introducer needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- One: ..... Introducer Needle: 21 Ga. x 2-1/4" (5.7 cm) TW with Echogenic Tip

**Rx only.**

**Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use — do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.



Arrow International, Inc.  
2400 Bernville Road  
Reading, PA 19605 USA


22 Ga. x 1-3/4" (4.45 cm)  
Catheter Assembly

.016" (.40 mm) dia. x  
13-1/8" (33 cm)  
Spring-Wire Guide

17 Ga. x 2-3/4"  
(7 cm) Peel-Away  
Sheath over 3 Fr. Dilator

21 Ga. x 1-1/2" (3.81 cm)  
TW Introducer Needle

21 Ga. x 2-1/4" (5.7 cm)  
TW Introducer Needle  
with Echogenic Tip



IF01351-103C (4/16)

# ARROW®

INTERNATIONAL

REF **PI-01351-LS5** to be used with Arrow **3** Fr. PICC **.016** inch dia. spring-wire guide  
**PICC Seldinger Conversion Set**

Peel to Open  
LATEX-FREE