

REF PI-01451-SS

Latex-Free

PICC Seldinger Conversion Set

for use with Arrow® 4 Fr. PICC

ARROW®

Contents:

- One: Sheath: 15 Ga. x 1-1/4" (3.18 cm) Radiopaque Peel-Away over 4 Fr. Dilator
- One: Spring-Wire Guide: .025" (.64 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends)
- One: Catheter: 18 Ga. x 1-9/16" (3.97 cm) Radiopaque over 20 Ga. RW Introducer Needle
- One: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.



Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA




18 Ga. x 1-9/16" (3.97 cm)
Catheter Assembly

15 Ga. x 1-1/4" (3.18 cm)
Peel-Away Sheath
over 4 Fr. Dilator

.025" (.64 mm) dia.
x 13-1/8" (33 cm)
Spring-Wire Guide

20 Ga. x 1-1/2"
(3.81 cm) XTW
Introducer Needle



1101451-101C (4/16)

REF
PI-01451-SS

for use with Arrow 4 Fr. PICC

.025 inch dia. spring-wire guide

PICC Seldinger Conversion Set

