## REF PI-01451-LS Lutex-Free PICC Seldinger Conversion Set for use with Arrow® 4 Fr. PICC ARRON(®

22 Ga. x 1-3/4" (4.45 cm) Catheter Assembly \_\_\_\_\_

**PICC Seldinger Conversion Set** 

Rx only.

Contents:

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

One:.... Introducer Needle: 21 Ga. x

One:..... Sheath: 15 Ga. x 2-3/4"

One:.....Spring-Wire Guide: .018" (.46 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip on Both Ends) One:....Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle One:....Introducer Needle: 21 Ga. x

(7 cm) Radiopaque Peel-Away over 4 Fr. Dilator

1-1/2" (3.81 cm) TW

2-3/4" (7 cm) TW with Echogenic Tip


Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains NO natural rubber latex.



Arrow International, Inc. 2400 Bernville Road, Reading, PA 19605 USA





Peel to Open UNITH-FREE