Rx only.

injury or death.

rubber.

Warning: Prior to use read all

package insert warnings, pre-

cautions, and instructions. Failure to do so may result in severe patient

Sterile. Do not use if package has

been previously opened or damaged.

Fluid path components are non-

pyrogenic. Contains no medication.

Single use — do not resterilize.

This product contains dry natural

Sterilized by ethylene oxide.

Arrow International, Inc. 2400 Bernville Road Reading, PA 19605 USA

REF AJ-09801

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port for use with 7.5 Fr. Catheters

Special Order Product

One: Sheath: 8.5 Fr. x 2-3/4" (7 cm) Arrow-Flex® Radiopaque Polyurethane

with Integral Hemostasis Valve/Side Port, Tissue Dilator and bonded High-Flow 3-way Stopcock

Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™1

Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)

One: Injection Needle: 25 Ga. x 1" (2.54 cm) One: Syringe: 3 mL Luer-Lock

Two: Syringes: 5 mL Luer-Slip

One: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter

One: Drape: 21" x 36" with 2-3/4" fenestration, with adhesive

Five: Gauze Pads: 4" x 4" One: Obturator: 8 Fr. One: Scalpel: #11 Two: Prep Sponge Swabs

One: Suture: 3-0 Silk with Straight Needle

One: Tray: Prep

One: SharpsAway® Disposal Cup

¹ US Patent No. 5,484,419.

