

REF AJ-09886

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port
for use with 7.5 - 8 Fr. Catheters

ARROW®

Special Order Product

Contents:

- One: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- One: Spring-Wire Guide, Marked: .025" (.64 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
- One: Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- One: Injection Needle: 25 Ga. x 1" (2.54 cm)
- One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
- One: Pressure Transduction Probe
- One: Syringe: 3 mL Luer-Lock
- One: Syringe: 5 mL Luer-Slip
- One: Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- One: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- One: Obturator: 8 Fr.
- One: SharpsAway® Disposal Cup
- One: CSR Wrap
- One: Drape: 24" x 36"
- One: Drape: 24" x 36" with 4" fenestration
- One: Scalpel: #11
- One: Stopcock: 3-way High-Flow
- Two: Gauze Pads: 2" x 2"
- Five: Gauze Pads: 4" x 4"
- Two: Prep Sponge Swabs
- One: Suture: 3-0 Braided Silk with Straight Needle
- One: Tray: Prep

¹ US Patent No. 5,484,419.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use – do not resterilize. Sterilized by ethylene oxide.

This product contains dry natural rubber.

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

U-09886-102B (9/10)

ARROW
INTERNATIONAL
Special Order Product

REF
AJ-09886

8.5 Fr.

10 cm
sheath
length

.025 inch dia.
spring-wire
guide

Percutaneous Sheath Introducer Kit

Peel to
Open