


**SHARPS
SAFETY
FEATURES**

PICC Kit with Placement Wire

REF **PK-01351-TWSSP**



INFECTION PROTECTION FEATURE:

- **ChloroPrep®** One-Step skin antiseptic, 2% chlorhexidine gluconate in 70% isopropyl alcohol with Hi-Lite Orange™ Tint

SAFETY FEATURES:

- **SharpsAway II™** Locking Disposal Cup
- **STATLOCK®** Catheter Stabilization Device
- **SafetyGlide™** Protected Needles

Catheter Diameter: 3 Fr., 0.99 mm

Priming Volume: 0.30 mL

Priming volumes are approximate and are done without the injection site caps.

Flow Rate: 150 mL/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Pump Flow Rate: 1040 mL/hr

Pump flow rates determined at maximum pump pressure of 10 psig.

Contents:

- 1: Indwelling Catheter, Trimmable: 3 Fr. x 55 cm Radiopaque Polyurethane with T-Port Connector¹, Integral Extension Line, Extension Line Clamps, Integral Suture Wing, Contamination Guard and Placement Wire**
- 1: GlideThru™ Peel-Away Sheath: 3 Fr. x 2-3/4" (7 cm) Radiopaque over 3 Fr. Dilator
- 1: Spring Wire Guide Assembly: .016" (0.40 mm) dia. x 13-1/8" (33 cm) (Straight Soft Tip On Both Ends)
- 1: Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle
- 1: Injection Needle: SafetyGlide™ 21 Ga. x 1-1/2" (3.81 cm) TW**
- 1: Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Lock
- 1: Syringe: 10 mL Luer-Slip
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 5/8" (1.60 cm) and 3 mL Luer-Lock Syringe**
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®² One-Step Solution with Hi-Lite Orange™ Tint**
- 1: Skin Protectant Prep Pad
- 1: SharpsAway II™ Locking Disposal Cup**
- 1: Catheter Trimmer
- 1: CSR Wrap
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Towel
- 1: Drape: 24" x 36"
- 1: Forceps
- 1: Dressing: OpSite®³ 14 cm x 10 cm
- 1: Dressing: STATLOCK®⁴ Catheter Stabilization Device**
- 5: Gauze Pads: 2" x 2"
- 6: Gauze Pads: 4" x 4"
- 2: Paper Tape Measures
- 2: Surgical Apparel: Masks
- 1: Tourniquet
- 1: Tape: Steri-Strip®⁵
- 1: Chart Sticker
- 1: Patient ID Card
- 1: Patient Information Booklet

¹ A trademark of Becton, Dickinson and Company.

² A registered trademark of CareFusion or one of its subsidiaries.

³ A registered trademark of T.J. Smith & Nephew, Limited.

⁴ A registered trademark of C. R. Bard, Inc.

⁵ A registered trademark of 3M Company.

Rx only.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterile. Single use – do not resterilize. Sterilized by ethylene oxide. Keep away from sunlight. Keep dry. Do not use if package has been previously opened or damaged. Not made with natural rubber latex. Fluid path components are non-pyrogenic. Contains medication.

***Precaution: Contains Phthalate: DEHP.**

Store between 15 - 30°C (59 - 86°F).

Avoid freezing and excessive heat above 40°C (104°F).



ARROW
INTERNATIONAL
2400 Bernville Road
Reading, PA 19605 USA
www.arrowvascular.com

IK-01351-105F (5/16)

**SHARPS
SAFETY
FEATURES**

PICC Kit with Placement Wire

3 Fr. / **55** cm catheter length / **.016** inch dia. spring-wire guide

REF **PK-01351-TWSSP**

Peel to
OPEN