



PICC Peripherally Inserted Central Catheter Kit with Blue FlexTip® Catheter and Placement Wire

REF **PK-01451-W**

4 Fr. | **55** cm catheter length

Contents:

- | | |
|---|---|
| <ul style="list-style-type: none"> 1: Indwelling Catheter: 4 Fr. x 55 cm Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Integral Extension Line, Extension Line Clamp, Integral Suture Wing, Contamination Guard and Placement Wire 1: Peel-Away Sheath: 15 Ga. x 1-1/4" (3.18 cm) Radiopaque over 17 Ga. TW Introducer Needle 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm) 1: Syringe: 10 mL Luer-Slip 1: Skin Protectant Prep Pad (non-sterile solution) 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution) 1: Dust Cap: Non-Vented 1: SecondSite™ Adjustable Hub: Catheter Clamp 1: SecondSite™ Adjustable Hub: Fastener 1: SharpsAway® Disposal Cup 1: Towel 1: Towel: 24" x 36" (61 cm x 91 cm) | <ul style="list-style-type: none"> 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration 1: Forceps 1: Patient ID Card 1: Chart Sticker 1: Patient Information Booklet 1: Dressing: OpSite®² 14 cm x 10 cm 2: Paper Tape Measure 1: Tourniquet 5: Gauze Pad: 2" x 2" (5 cm x 5 cm) 6: Gauze Pad: 4" x 4" (10 cm x 10 cm) 1: Dressing: STATLOCK®³ Catheter Stabilization Device 2: Surgical Apparel: Mask 1: Surgical Tape: Transpore™ 1" x 54" Roll |
|---|---|

¹A registered trademark of CareFusion or one of its subsidiaries.
²A registered trademark of T.J. Smith & Nephew, Limited.
³A registered trademark of C. R. Bard, Inc.

DuPont™ Tyvek®
polypropylene

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.

Contraindications: This device is contraindicated wherever there is presence of device related infections, previous or current thrombosis. Clinical assessment of patient must be completed to ensure no contraindications exists.



STERILE EO




~~LATEX~~

Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (18 Ga.)	0.52	850	3700

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig.

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IK-01451-1021 (7/17)



PICC Peripherally Inserted Central Catheter Kit with Blue FlexTip® Catheter and Placement Wire

REF **PK-01451-W**

4 Fr. | **55** cm catheter length