



Two-Lumen Peripherally Inserted Central Catheter Kit with Blue FlexTip® Catheter and Placement Wire

REF **PK-01452-W**

4 Fr. | **2** Lumen | **55** cm catheter length

Contents:

- 1: Two-Lumen Indwelling Catheter: 4 Fr. x 55 cm Radiopaque Polyurethane with Blue FlexTip®, T-Port Connector, Extension Line Clamps, Contamination Guard and Placement Wire
- 1: Peel-Away Sheath: 15 Ga. x 1-1/4" (3.18 cm) Radiopaque over 17 Ga. TW Introducer Needle
- 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Syringe: 10 mL Luer-Slip
- 1: Skin Protectant Prep Pad (non-sterile solution)
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 2: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SharpsAway® Disposal Cup
- 1: Towel
- 1: Towel: 24" x 36" (61 cm x 91 cm)
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Forceps
- 1: Patient ID Card
- 1: Chart Sticker
- 1: Patient Information Booklet
- 1: Dressing: OpSite®² 14 cm x 10 cm
- 2: Paper Tape Measure
- 1: Tourniquet

- 5: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 6: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Dressing: STATLOCK®³ Catheter Stabilization Device
- 2: Surgical Apparel: Mask
- 1: Surgical Tape: Transpore™⁴ 1" x 54" Roll

¹A registered trademark of CareFusion or one of its subsidiaries.
²A registered trademark of T.J. Smith & Nephew, Limited.
³A registered trademark of C. R. Bard, Inc.
⁴A registered trademark of 3M Company.


Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



Lumen	Priming Volume* (mL)	Flow Rate† (mL/hr)	Pump Flow Rate‡ (mL/hr)
Distal (22 Ga.)	0.30	125	680
Proximal (22 Ga.)	0.40	150	900

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.
 ‡ Pump flow rates are determined at maximum pump pressure of 10 psig.

Not made with natural rubber latex

Store between 15 - 30°C (59 - 86°F).
 Avoid freezing and excessive heat above 40°C (104°F)

Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA
Teleflex®
www.teleflex.com

IK-01452-101J (07/19)



Two-Lumen Peripherally Inserted Central Catheter Kit with Blue FlexTip® Catheter and Placement Wire

REF **PK-01452-W**

4 Fr. | **2** Lumen | **55** cm catheter length **ARROW®**