

Product No. AK-04300

# Central Venous Catheterization Kit

# ARROW®

Now Contains  
Tinted  
ChloraPrep®  
Hi-Lite  
Orange™

Catheter Diameter: 16 Ga., 1.7 mm

Priming Volume: 0.22 cc

Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 cc.

Flow Rate: 3400 cc/hr

Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

### Contents:

- One: . . . . . Indwelling Catheter: 16 Ga. x 8" (20 cm)  
Radiopaque Polyurethane with Integral Suture Wing
- One: . . . . . Spring-Wire Guide: .032" (.81 mm) dia. x  
17-13/16" (45 cm) (Straight Soft Tip on  
One End - "J" Tip on Other)
- One: . . . . . Catheter: 18 Ga. x 2-1/2" (6.35 cm)  
Radiopaque over 20 Ga. RW Introducer Needle
- One: . . . . . Injection Needle: 20 Ga. x 1-1/2" (3.81 cm)  
and 5 cc Luer-Slip Syringe
- One: . . . . . Injection Needle: 25 Ga. x 1" (2.54 cm)  
and 3 cc Luer-Lock Syringe
- One: . . . . . Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW  
and 5 cc Luer-Slip Syringe
- One: . . . . . Tissue Dilator
- One: . . . . . 3 mL Applicator Pouch, 2% CHG and  
70% IPA ChloraPrep®<sup>1</sup> One-Step Solution  
with Hi-Lite Orange™ Tint
- One: . . . . . .5 mL HCl, 1% Lidocaine Solution
- One: . . . . . Clamp: Catheter
- One: . . . . . Fastener: Catheter Clamp
- One: . . . . . SharpsAway® Disposal Cup
- Three: . . . . . Drapes: 26" x 18" with adhesive
- One: . . . . . Local Anesthetic Components
- One: . . . . . Scalpel: #11
- Three: . . . . . Gauze Pads: 4" x 4"
- One: . . . . . Suture: 3-0 Nylon with Straight Needle
- One: . . . . . Tubing: 6" Extension±

<sup>1</sup> A registered trademark of Enturia, Inc.

### Rx only.

**Warning:** Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated. Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

<sup>±</sup>Precaution: Contains Phthalate: DEHP.

Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Not made with natural rubber latex.



Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA

11-04300-100K (3/16)

Product No. **AK-04300** **16**<sup>Ga.</sup> **20** cm catheter length **.032** inch dia. spring-wire guide  
Central Venous Catheterization Kit



Product No. **AK-04300** **16**<sup>Ga.</sup> **20** cm catheter length **.032** inch dia. spring-wire guide  
Central Venous Catheterization Kit

Peel to Open