





Arterial Line Kit

20^{Ga.} 12^{cm} catheter length .02

- Indwelling Catheter: 20 Ga. x 12 cm Radiopaque FEP Spring-Wire Guide: .025" (0.64 mm) dia. x 13-3/4" (35 cm) (Straight Soft Tip on One End "J" Tip on Other) Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle

- 1: Injection Needle: 25 Ga. x 1" (2.54 cm) 1: Injection Needle: 18 Ga. x 1-1/2" (3.81 cm) 1: Introducer Needle: 20 Ga. x 2" (5.08 cm) XTW 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

- 2: Syringe: 3 mL Luer-Lock
 1: 5 mL 1% Lidocaine HCI Solution and Alcohol Prep
 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
 Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- 3: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

⚠ California Prop. 65
Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.





















Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-04510-105A (2020-02)







