

REF AK-05801

ARROW®

Trauma Kit

8.5^{Fr.} | 8.89^{cm} sheath length | .035^{inch dia.} spring-wire guide

Contents:

- 1: Sheath: 8.5 Fr. x 8.89 cm Arrow-Flex® Radiopaque FEP with Tissue Dilator
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Syringe: 5 mL Luer-Slip
- 1: Syringe: 30 mL Luer-Lock
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloroPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 1: Dressing: OpSite®² 14 cm x 10 cm
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

²A registered trademark of T.J. Smith & Nephew, Limited.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.



STERILE EO



LATEX
Not made with natural rubber latex

20°C (68°F) | 25°C (77°F)
Store between 20-25°C (68-77°F).
Avoid excessive heat above 40°C (104°F).

REF AK-05801

ARROW®

Trauma Kit

8.5^{Fr.} | 8.89^{cm} sheath length | .035^{inch dia.} spring-wire guide

Manufacturer:
Arrow International, Inc.
Subsidiary of Teleflex Incorporated
2400 Bernville Road
Reading, PA 19605 USA

Teleflex®
www.teleflex.com

IK-05801-105A (2020/04)