

REF AK-09000

ARROW®

# Peritoneal Lavage Kit

8<sup>Fr.</sup> | 24<sup>cm</sup> catheter length | .035<sup>inch dia.</sup> spring-wire guide

## Contents:

- 1: Peritoneal Lavage Catheter: 8 Fr. x 24 cm Radiopaque FEP with Side Holes
- 1: Spring-Wire Guide: .035" (0.89 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Injection Needle: SafetyGlide™ 25 Ga. x 1" (2.54 cm)
- 1: Injection Needle: SafetyGlide™ 22 Ga. x 1-1/2" (3.81 cm) RW
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 2: Syringe: 5 mL Luer-Slip
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint
- 1: SharpsAway II™ Locking Disposal Cup
- 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Safety Scalpel: #11
- 1: Medication Label: 1% Lidocaine
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 4: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Tubing Set: Fluid Administration

<sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

⚠ **California Prop. 65**

**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.

**Warning: The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male premature infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.**

Arrow International provides the enclosed medication label(s) for your convenience. Please ensure that these labels are applied to the correct syringe and corresponding medication.



~~LATEX~~  
Not made with natural rubber latex

PHT DEHP  
Precaution:  
Contains Phthalate: DEHP

20°C (68°F) - 25°C (77°F)  
Store between 20-25°C (68-77°F).

REF AK-09000

ARROW®

# Peritoneal Lavage Kit

8<sup>Fr.</sup> | 24<sup>cm</sup> catheter length | .035<sup>inch dia.</sup> spring-wire guide

Manufacturer:  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)

IK-09000-105A, Rev. 01 (2020-12)