

Antimicrobial surfaces utilizing
chlorhexidine acetate and
silver sulfadiazine

Contents:

- One: Sheath: 8.5 Fr. x 4" (10 cm) Arrow-Flex®
Radiopaque Polyurethane with ARROWgard®
Antimicrobial Surface Treatment¹, Integral Hemostasis
Valve/Side Port and Tissue Dilator
- One: Spring-Wire Guide: .035" (.89 mm) dia. x
17-13/16" (45 cm) (Straight Soft Tip
on One End – "J" Tip on Other)
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm)
Radiopaque over 20 Ga. RW Introducer Needle
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
and 5 cc Luer-Slip Syringe
- One: Introducer Needle: 18 Ga. x
2-1/2" (6.35 cm) XTW and 5 cc Luer-Slip Syringe
- One: Injection Needle: 25 Ga. x 1" (2.54 cm)
and 3 cc Luer-Lock Syringe
- One: Suture: 3-0 Braided Silk with Straight Needle
- One: 5 ml HCl, 1% Lidocaine Solution
- One: 1 Applicator Pouch, 2% CHG and
70% IPA Chloraprep®² One-Step Solution
- One: Dust Cap: Non-Vented
- One: Cath-Gard® Catheter Contamination Shield:
80 cm with Sterile Tape
- One: Drape: 24" x 36"
- One: Drape: 24" x 36" with 4" fenestration
- One: Local Anesthetic Components
- One: Dressing: Tegaderm®³ 10 cm x 12 cm
- Two: Gauze Pads: 2" x 2"
- Five: Gauze Pads: 4" x 4"
- One: Obturator: 8 Fr.
- One: Scalpel: #11
- One: Stopcock: 4-way High-Flow
- One: Tubing: 20" Extension=
- One: SharpsAway® Disposal Cup

¹Licensed under US Patent Nos. 4,581,028 and 5,019,096.

²A registered trademark of Medi-Flex Hospital Products Corporation.

³A registered trademark of 3M Company.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

Contraindications: The ARROWgard Blue® antimicrobial sheath introducer is contraindicated for patients with known hypersensitivity to chlorhexidine acetate, silver sulfadiazine, and/or sulfa drugs (refer to product instructions for references).

All contents sterile unless otherwise indicated.

Single use – do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

±Precaution: Contains Phthalate: DEHP.

*Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Not made with natural rubber latex.

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

**INFECTION
PROTECTION
FEATURE:**

- **Chloraprep®***
One-Step skin antiseptic,
2% chlorhexidine
gluconate in 70%
isopropyl alcohol



Product No. AK-09804-AG

Percutaneous Sheath Introducer Kit

with ARROWgard Blue® Sheath
and Integral Hemostasis Valve/Side Port
for use with 7 - 7.5 Fr. Catheters

ARROW®

Product No. **AK-09804-AG** **8.5** Fr. **10** cm sheath length **.035** inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with ARROWgard Blue® Sheath

Product No. **AK-09804-AG** **8.5** Fr. **10** cm sheath length **.035** inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit with ARROWgard Blue® Sheath

IK09804-103, Rev.1 (10/16)

Peel to Open