

Product No. AK-09903-S

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port
for use with 7.5 - 8 Fr. Catheters

ARROW®

Contents:

- One: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
 - One: Spring-Wire Guide, Marked: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End – "J" Tip on Other) with Arrow Advancer™¹
 - One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
 - One: Pressure Transduction Probe
 - One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm)XTW and 5 cc Arrow® Raulerson Spring-Wire Introduction Syringe²
 - One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 cc Luer-Slip Syringe
 - One: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 cc Luer-Lock Syringe
 - One: Suture: 3-0 Braided Silk with Straight Needle
 - One: .5 ml HCl, 1% Lidocaine Solution
 - One: 1 Applicator Pouch, 2% CHG and 70% IPA ChloraPrep®³ One-Step Solution
 - One: Dust Cap: Non-Vented
 - One: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
 - One: Drape: 24" x 36" with 4" fenestration
 - One: Drape: 24" x 36"
 - One: Local Anesthetic Components
 - Two: Gauze Pads: 2" x 2"
 - Five: Gauze Pads: 4" x 4"
 - One: Obturator: 8 Fr.
 - One: Scalpel: #11
 - One: Stopcock: 4-way High-Flow
 - One: Tubing: 20" Extension⁴
 - One: SharpsAway® Disposal Cup
- ¹US Patent Nos. 5,484,419 and 6,477,402.
²Licensed under US Patent Nos. 4,813,938 and 5,045,065.
³A registered trademark of Medi-Flex Hospital Products Corporation.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use — do not resterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

Not made with natural rubber latex.

*Precaution: Contains Phthalate: DEHP.

*Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605 USA

**INFECTION
PROTECTION
FEATURE:**

- ChloraPrep®*
One-Step skin antiseptic,
2% chlorhexidine
gluconate in 70%
isopropyl alcohol

ARROW®

INTERNATIONAL

Product No.
AK-09903-S

9^{Fr.}

10^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

Percutaneous Sheath Introducer Kit

Product No.
AK-09903-S

9^{Fr.}

10^{cm}
sheath
length

.035^{inch dia.}
spring-wire
guide

Percutaneous Sheath Introducer Kit

IK09903-1000, Rev. 1 (10/16)

Peel to
Open