

REF AK-09903-A

Latex-Free

# Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port  
for use with 7.5 - 8 Fr. Catheters

# ARROW®

**Contents:**

- One: Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex® Radiopaque Polyurethane with Integral Hemostasis Valve/Side Port and Tissue Dilator
- One: Spring-Wire Guide, Marked: .035" (.89 mm) dia. x 17-13/16" (45 cm) (Straight Soft Tip on One End — "J" Tip on Other) with Arrow Advancer™
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- One: Pressure Transduction Probe
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- One: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe
- One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe<sup>2</sup>
- One: 3 mL Applicator Pouch, 2% CHG and 70% IPA Chloraprep® One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- One: Cath-Gard® Catheter Contamination Shield: 80 cm with TwistLock™ Adapter
- One: Obturator: 8 Fr.
- One: Suture: 3-0 Braided Silk with Straight Needle
- One: Dust Cap: Non-Vented
- One: SharpsAway® Disposal Cup
- One: Drape: 24" x 36"
- One: Drape: 24" x 36" with 4" fenestration
- One: CSR Wrap
- One: Scalpel: #11
- One: Stopcock: 4-way High-Flow
- Seven: Gauze Pads: 4" x 4"
- One: Tubing: 20" Extension

<sup>1</sup>US Patent No. 5,484,419.  
<sup>2</sup>Licensed under US Patent No. 5,045,065.  
<sup>3</sup>A registered trademark of Cardinal Health, Inc., or one of its subsidiaries.

**Rx only.**

**Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.**

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication. Single use — do not resterilize. Sterilized by ethylene oxide.

\*Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

This product contains **NO** natural rubber latex.


Arrow International, Inc.  
2400 Bernville Road, Reading, PA 19605 USA

www.teleflex.com



**INFECTION PROTECTION FEATURE:**

- Chloraprep®\*  
One-Step skin antiseptic,  
2% chlorhexidine gluconate in  
70% isopropyl alcohol



IK09903-1158 (11/17)

**ARROW®**  
INTERNATIONAL

REF  
AK-09903-A

9 Fr.

10 cm sheath length

.035 inch dia. spring-wire guide

Percutaneous Sheath Introducer Kit

Peel to Open  
LATEX-FREE