

Product No. AK-09907-S

Percutaneous Sheath Introducer Kit

with Integral Hemostasis Valve/Side Port
for use with 7.5 Fr. - 8 Fr. Catheters

ARROW®

Contents:

- One:..... Sheath: 9 Fr. x 4" (10 cm) Arrow-Flex®
Radiopaque Polyurethane with
Integral Hemostasis Valve/Side Port and Tissue Dilator
 - One:..... Spring-Wire Guide, Marked: .025" (.64 mm) dia.
x 17-13/16" (45 cm)
(Straight Soft Tip on One End – "J" Tip
on Other) with Arrow Advancer™¹
 - One:..... Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque
over 22 Ga. TW Introducer Needle
 - One:..... Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) RW
 - One:..... Pressure Transduction Probe
 - One:..... Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW
and 5 cc Arrow® Raulerson Spring-Wire
Introduction Syringe²
 - One:..... Injection Needle: 22 Ga. x 1-1/2" (3.81 cm)
and 5 cc Luer-Slip Syringe
 - One:..... Injection Needle: 25 Ga. x 1" (2.54 cm)
and 3 cc Luer-Lock Syringe
 - One:..... 5 ml HCl, 1% Lidocaine Solution
 - One:..... 1 Applicator Pouch, 2% CHG and
70% IPA ChloraPrep®³ One-Step Solution
 - One:..... Dust Cap: Non-Vented
 - One:..... Cath-Gard® Catheter Contamination Shield:
80 cm with TwistLock™ Adapter
 - One:..... Drape: 24" x 36" with 4" fenestration
 - One:..... Drape: 24" x 36"
 - One:..... Local Anesthetic Components
 - Two:..... Gauze Pads: 2" x 2"
 - Five:..... Gauze Pads: 4" x 4"
 - One:..... Obturator: 8 Fr.
 - One:..... Scalpel: #11
 - One:..... Stopcock: 4-way High-Flow
 - One:..... Suture: 3-0 Silk with Straight Needle
 - One:..... Tubing: 20" Extension±
 - One:..... SharpsAway® Disposal Cup
- ¹US Patent Nos. 5,484,419 and 6,477,402.
²Licensed under US Patent Nos. 4,813,938 and 5,045,065.
³A registered trademark of Medi-Flex Hospital Products Corporation.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Single use – do not re-sterilize. Sterilized by ethylene oxide. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains medication.

±Precaution: Contains Phthalate: DEHP

*Store between 20 - 25°C (68 - 77°F). Avoid freezing and excessive heat above 40°C (104°F).

Not made with natural rubber latex.

Arrow International, Inc.
2400 Bernville Road
Reading, PA 19605

**INFECTION
PROTECTION
FEATURE:**

- **ChloraPrep®***
One-Step skin antiseptic,
2% chlorhexidine
gluconate in 70%
isopropyl alcohol

ARROW
INTERNATIONAL

Product No. **AK-09907-S** **9** Fr. **10** cm sheath length **.025** inch dia. spring-wire guide
Percutaneous Sheath Introducer Kit

Product No. **AK-09907-S** **9** Fr. **10** cm sheath length **.025** inch dia. spring-wire guide
Percutaneous Sheath Introducer Kit

IK-09907-101H, Rev. 1 (10/16)

Peel to
Open