



Multi-Lumen Hemodialysis Catheterization Kit  
for High Volume Infusions

REF **AK-12123-PCDC**      **12** Fr. | **3** Lumen | **16** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 16 cm Pressure Injectable Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps

1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer

1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle

1: Injection Needle: SafetyGlide™<sup>1</sup> 25 Ga. x 1" (2.54 cm)

1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

1: Pressure Transduction Probe

1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe

1: Injection Needle: SafetyGlide™<sup>1</sup> 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe

1: Syringe: 3 mL Luer-Lock

1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm

1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)

1: SharpsAway II™ Locking Disposal Cup

1: SharpsAway® Disposal Cup

1: Drape: 41" x 55" with 4" fenestration

1: Particulate Filter: 5 Micron Straw

1: Needle Holder

1: Safety Scalpel: #11

1: Flow Rate Information Card with Injection Log

2: Gauze Pad: 2" x 2" (5 cm x 5 cm)

5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

1: Suture: 3-0 Silk with Curved Needle
- <sup>1</sup>A trademark of Becton, Dickinson and Company.

<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.**  
Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
⚠️ **California Prop. 65**  
**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Fluid path components are non-pyrogenic.  
**Contraindications: The Multi-Lumen Acute Hemodialysis Catheter (AHDC) is not designed for long-term (≥ 30 days) hemodialysis or for insertion into thrombosed vessels.**



🏢 **Manufacturer:**  
Arrow International, Inc.  
Subsidiary of Teleflex Incorporated  
2400 Bernville Road  
Reading, PA 19605 USA  
**Teleflex®**  
[www.teleflex.com](http://www.teleflex.com)  
IK-12123-103C (9/18)



Multi-Lumen Hemodialysis Catheterization Kit  
for High Volume Infusions

REF **AK-12123-PCDC**      **12** Fr. | **3** Lumen | **16** cm catheter length | **.035** inch dia. spring-wire guide      **ARROW®**