



Multi-Lumen CVC Kit with Blue FlexTip® Catheter

REF AK-12703-CDC

7 Fr. 3 Lumen 16 cm catheter length .032

Contents: _

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer™
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- **Pressure Transduction Probe**
- Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 mL Luer-Lock Syringe 1:
- Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Local Anesthetic Components 1:
- 1: Tissue Dilator
- 5 mL 1% Lidocaine HCl Solution


- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®1 One-Step Solution with Hi-Lite Orange[™] Tint
- Dust Cap: Non-Vented
- SecondSite[™] Adjustable Hub: Fastener 1:
- SecondSite[™] Adjustable Hub: Catheter Clamp
- SharpsAway II™ Locking Disposal Cup
- Drape: 41" x 55" with 4" fenestration 1:
- 1: Safety Scalpel: #11
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.



Caution



instructions for use







Sterilized by ethylene oxide








package is damaged



natural rubber latex



Store between 20 - 25°C (68 - 77°F).



Lumen	Priming Volume* (mL)	Gravity Flow Rate [†] (mL/hr.)
Distal (16 Ga.)	0.40	3550
Medial (18 Ga.)	0.40	1870
Proximal (18 Ga.)	0.40	2050

Manufacturer: Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-12703-107F (1/17)

ARRO



Multi-Lumen CVC Kit with Blue FlexTip® Catheter



Priming volumes are approximate and are done without accessories.
Flow rates are done with deionized water or norma saline, room temperature, 100 cm head height and represent approximate flow capabilities.