





Pediatric Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

2^{Lumen} 3 cm catheter length **5**^{Fr.} inch dia. spring-wire guide

Contents

- 1: Two-Lumen Indwelling Catheter: 5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Spring-Wire Guide: .021" (0.53 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: 1: 1:
- Injection Needle: 25 Ga. x 1-7(2:34 cm) Introducer Needle: 20 Ga. x 1-1/2" (3.81 cm) XTW Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm) Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Slip 1:
- 1:
- 1:
- 1:
- Catheter: 20 Ga. x 1-3/4" (4.45 cm) Radiopaque over 22 Ga. TW Introducer Needle with 5 mL Luer-Syringe Syringe: 3 mL Luer-Lock Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm 5 mL 1% Lidocaine HCI Solution and Alcohol Prep 3 mL Applicator 2% CHG and 70% IPA ChloraPrep^{®1} One-Step Solution with Hi-Lite Orange[™] Tint (non-sterile solution) Dust Can: Non-Vented 1: 1:
- 2: Dust Cap: Non-Vented

- 2: Dust Cap: Non-Vented
 1: SecondSite[™] Adjustable Hub: Fastener
 1: SecondSite[™] Adjustable Hub: Catheter Clamp
 1: SharpsAway[®] Disposal Cup
 1: Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
 1: Scalpel: #11
 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
 3: Gauze Dad: 2" x 4" (10 cm x 10 cm)

- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Sterilized using ethylene oxide unless otherwise indicated in the contents list.

2

STERINKE

\land California Prop. 65

 \wedge

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov

Fluid path components are non-pyrogenic.

 $\mathbf{\tilde{i}}$


1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.



*

STERILE EO



Priming volumes are approximate and are done without accessories. Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

 $(\bigotimes$



LÀTEX

Not made with natural rubber latex



Manufacturer Arrow International, Inc Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-14502-135A (2020-02)

REF AK-14502



Pediatric Two-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter 5^{Fr.} 2^{Lumen} 13^{cm}_{catheter} .021^{inch dia.} spring-wire guide

