



## Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter For Femoral Vein Insertion

**REF AK-14553** 

## **Contents:**

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 30 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- Spring-Wire Guide: .018" (0.46 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- Injection Needle: 25 Ga. x 1" (2.54 cm)
- Introducer Needle: 21 Ga. x 1-1/2" (3.81 cm) TW
- Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- Catheter: 22 Ga. x 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Syringe: 3 mL Luer-Lock
- Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm
- Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm
- 5 mL 1% Lidocaine HCl Solution and Alcohol Prep (non-sterile solution)
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®1 One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint (non-sterile solution)

- Dust Cap: Non-Vented
- SecondSite™ Adjustable Hub: Fastener SecondSite™ Adjustable Hub: Catheter Clamp
- SharpsAway® Disposal Cup
- Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- Scalpel: #11
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
- Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Suture: 3-0 Braided Silk with Straight Needle

<sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. Warning: Read all package insert warnings, cautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

**↑** California Prop. 65

Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Fluid path components are non-pyrogenic.




















Not made with natural rubber latex



Store between 20 - 25°C (68 - 77°F).



Lumen	Priming Volume* (mL)	Flow Rate <sup>†</sup> (mL/hr.)
Distal (20 Ga.)	0.39	370
Medial (22 Ga.)	0.38	320
Proximal (22 Ga.)	0.41	320

Priming volumes are approximate and are done

Arrow International, Inc. Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-14553-130D (01/19)



Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter For Femoral Vein Insertion



<sup>†</sup> Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.