

Product No. AK-14703-J

LATEX-FREE


Arrow-Howes™ Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

ARROW®

Lumen	Priming Volume* (cc)	Flow Rate† (cc/hr)
Distal (16 Ga.)	0.49	2300
Medial (18 Ga.)	0.40	1000
Proximal (18 Ga.)	0.44	1100

*Priming volumes are approximate and are done without the injection cap. Injection cap priming volume is 0.17 cc.

† Flow rates are done with normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Contents:

- One: Multi-Lumen Indwelling Catheter: 7 Fr. x 12" (30 cm) Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps, Injection Site Caps
- One: Spring-Wire Guide, Marked: .032" (.81 mm) dia. x 26-13/16" (68 cm) (Straight Soft Tip on One End—"J" Tip on Other) with Arrow Advancer™¹
- One: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- One: Fastener: Catheter Clamp
- One: Pressure Transduction Probe
- One: Introducer Needle: 22 Ga. x 2-3/4" (6.98 cm) TW
- One: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 cc Arrow Raulerson Spring-Wire Introduction Syringe²
- One: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 cc Luer-Slip Syringe
- One: Injection Needle: 25 Ga. x 1" (2.54 cm) and 3 cc Luer-Lock Syringe
- Three: Dust Caps: Non-Vented
- One: Tissue Dilator
- One: CSR Wrap
- One: Drape: 21" x 36" with 2-3/4" fenestration, with adhesive
- Two: Gauze Pads: 2" x 2"
- Five: Gauze Pads: 4" x 4"
- One: Scalpel: #11
- Two: Prep Sponge Swabs
- Three: Stopcocks: 3-way
- One: Suture: 3-0 Silk with Straight Needle
- One: Tray: Prep
- Three: UserGard® Needle-Free Injection Hubs
- One: SharpsAway® Disposal Cup
- One: Clamp: Catheter

¹US Patent Nos. 5,484,419 and 6,477,402.


²Licensed under US Patent Nos. 4,813,938 and 5,045,065.

Rx only.

Warning: Prior to use read all package insert warnings, precautions, and instructions. Failure to do so may result in severe patient injury or death. Sterile. Do not use if package has been previously opened or damaged. Fluid path components are non-pyrogenic. Contains no medication. Single use — do not resterilize. Sterilized by ethylene oxide.

This product contains **NO** natural rubber latex.

Arrow International, Inc.
2400 Bernville Road, Reading, PA 19605 USA



IK14703-103L (5/03)



Product No. **AK-14703-J** **7 Fr.** **3 Lumen** **30 cm** catheter length **.032** inch dia. spring-wire guide

Arrow-Howes™ Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®

Peel to Open
LATEX-FREE