

REF **AK-15123-F**



Arrow-Howes™ Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter for High Volume Infusions

12 Fr. | **3** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 26-3/4" (68 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®¹ One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 3: Dust Cap: Non-Vented
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

⚠ **California Prop. 65**


Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

Sterilized using ethylene oxide unless otherwise indicated in the contents list.

Fluid path components are non-pyrogenic.

Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/min)
Distal (16 Ga.)	0.40	80
Medial (12 Ga.)	1.60	285
Proximal (12 Ga.)	1.50	305

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.



Not made with natural rubber latex



Store between 20-25°C (68-77°F).

REF **AK-15123-F**



Arrow-Howes™ Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter for High Volume Infusions

12 Fr. | **3** Lumen | **20** cm catheter length | **.035** inch dia. spring-wire guide

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA



www.teleflex.com

IK-15123-120A (2020-03)