



Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **AK-15703**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide

Contents:

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps
- 1: Spring-Wire Guide: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Syringe: 3 mL Luer-Lock
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle with 5 mL Luer-Slip Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Local Anesthetic Components
- 1: Tissue Dilator
- 1: 5 mL 1% Lidocaine HCl Solution
- 1: 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange™ Tint
- 3: Dust Cap: Non-Vented
- 1: SecondSite™ Adjustable Hub: Fastener
- 1: SecondSite™ Adjustable Hub: Catheter Clamp
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" with 4" fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle

¹A registered trademark of CareFusion or one of its subsidiaries.

DuPont™ Tyvek®
packaging

Rx only. Warning: Read all package insert warnings, precautions and instructions prior to use. Failure to do so may result in severe patient injury or death.

All contents sterile unless otherwise indicated.

Warning: This product contains a chemical known to the State of California to cause cancer, birth defects or other reproductive harm.

Fluid path components are non-pyrogenic.

Caution

Consult instructions for use

Do not reuse

Do not resterilize

Sterilized by ethylene oxide


Keep away from sunlight

Keep dry

Do not use if package is damaged

Not made with natural rubber latex

Store between 20 - 25°C (68 - 77°F)



Lumen	Priming Volume* (mL)	Gravity Flow Rate† (mL/hr.)
Distal (16 Ga.)	0.47	3200
Medial (18 Ga.)	0.40	1500
Proximal (18 Ga.)	0.42	1600

* Priming volumes are approximate and are done without accessories.
 † Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.

Manufacturer:
 Arrow International, Inc.
 Subsidiary of Teleflex Incorporated
 2400 Bernville Road
 Reading, PA 19605 USA

Teleflex®
 www.teleflex.com

IK-15703-1005 (2/17)



Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter

REF **AK-15703**

7 Fr. | **3** Lumen | **20** cm catheter length | **.032** inch dia. spring-wire guide