## REF AK-15703-ACDC



## Multi-Lumen CVC Kit with Blue FlexTip® Catheter



## Contents

- 1: Multi-Lumen Indwelling Catheter: 7 Fr. x 20 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps Spring-Wire Guide, Marked: .032" (0.81 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on
- 1: Other) with Arrow Advancer
- Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
   Injection Needle: 25 Ga. x 1" (2.54 cm)
   Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)

- 1.
- Pressure Transduction Probe Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction 1:
- Syringe 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- Syringe: 3 mL Luer-Lock Tissue Dilator: 8.5 Fr. (2.8 mm) x 10.2 cm 1: 1:
- 1.
- 3 mL 1% Lidocaine HCI Solution and Alcohol Prep 3 mL Applicator 2% CHG and 70% IPA ChloraPrep<sup>®1</sup> One-Step Solution with Hi-Lite Orange<sup>™</sup> Tint 1:
- (non-sterile solution)


- Inon-sterile solution)
  Dust Cap: Non-Vented
  SecondSite<sup>™</sup> Adjustable Hub: Fastener
  SecondSite<sup>™</sup> Adjustable Hub: Catheter Clamp
  SharpsAway II<sup>™</sup> Locking Disposal Cup
  SharpsAway<sup>®</sup> Disposal Cup
  Drape: 41" x 55" with 4" fenestration
  Safety Scalpel: #11

Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. Sterilized using ethylene oxide unless otherwise indicated in the contents list.


 $\wedge$ 

,----, , ,

⚠ California Prop. 65 Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov. Fluid path components are non-pyrogenic.



- Gauze Pad: 2 × 2 (Schrx 5 chr)
   Gauze Pad: 4" x 4" (10 cm x 10 cm)
   Suture: 3-0 Braided Silk with Straight Needle
- <sup>1</sup>A registered trademark of CareFusion or one of its subsidiaries.



Priming Volume\* Gravity Flow Rate<sup>†</sup> (mL) (mL/hr) 0.47 3200 Medial (18 Ga.) Proximal (18 Ga.) 0.40 0.42 1500 1600

ARROW

Priming volumes are approximate and are done without accessories. Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.









 $(\mathfrak{D})$ 



Manufacturer Arrow International, Inc Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-15703-130A (12/19)





Multi-Lumen CVC Kit with Blue FlexTip® Catheter

7<sup>Fr.</sup> 3<sup>Lumen</sup> 20<sup>cm</sup><sub>catheter</sub> 032<sup>inch dia.</sup>