





Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter



Contents

- 1: Multi-Lumen Indwelling Catheter: 5.5 Fr. x 13 cm Radiopaque Polyurethane with Blue FlexTip®, Extension Line Clamps Spring-Wire Guide: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - "J" Tip on Other)
- 1:
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)

- Inforction Needle: 25 Ga. X 1 (2:34 cm)
 Introducer Needle: 21 Ga. X 1-1/2" (3.81 cm) TW
 Blunt Fill Needle: 18 Ga. X 1-1/2" (3.81 cm)
 Catheter: 22 Ga. X 1-3/4" (4.45 cm) Radiopaque over 25 Ga. RW Introducer Needle with 5 mL Luer-Slip
- Syringe Injection Needle: 25 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe 1:
- Syringe: 3 mL Luer-Lock Tissue Dilator: 6 Fr. (2.0 mm) x 7.6 cm Tissue Dilator: 6 Fr. (2.0 mm) x 4.1 cm 5 mL 1% Lidocaine HCI Solution and Alcohol Prep 1:
- 1:
- 1: 1:
- 3 mL Applicator 2% CHG and 70% IPA ChloraPrep®¹ One-Step Solution with Hi-Lite Orange[™] Tint Dust Cap: Non-Vented SecondSite[™] Adjustable Hub: Fastener SecondSite[™] Adjustable Hub: Catheter Clamp 1.
- 1:
- 1:
- 1: 1:
- SharpsAway[®] Disposal Cup Drape: 18" x 26" (46 cm x 66 cm) with 3" (7.6 cm) fenestration
- 1: Scalpel: #11
- Gauze Pad: 2" x 2" (5 cm x 5 cm) Gauze Pad: 4" x 4" (10 cm x 10 cm)
- Gauze Pad: 2" x 2" (5 cm x 5 cm)
 Gauze Pad: 4" x 4" (10 cm x 10 cm)
 Suture: 3-0 Braided Silk with Straight Needle


Rx only. Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. ⚠ California Prop. 65 Warning: Cancer and Reproductive Harm www.P65Warnings.ca.gov.

2

STERINKE

i

Fluid path components are non-pyrogenic.



¹A registered trademark of CareFusion or one of its subsidiaries

| n | Priming Volume* (mL) | Gravity Flow Rate [†] (mL/hr) | |
|--|----------------------------|---|--|
| (20 Ga.) | 0.30 | 1200 | |
| l (22 Ga.) | 0.31 | 800 | |
| nal (22 Ga.) | 0.34 | 900 | |
| ning volumes are approximate and are done without acce | | | |

Priming volumes are approximate and are done without accessories.
 Flow rates are done with deionized water or normal saline, room temperature, 100 cm head height and represent approximate flow capabilities.









(&)



Manufacturer Arrow International, Inc Subsidiary of Teleflex Incorporated 2400 Bernville Road Reading, PA 19605 USA



IK-16553-131A (2020/04)

REF AK-16553



Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip® Catheter 5.5^{Fr.} 3^{Lumen} 13^{cm}_{catheter} 018^{inch dia.}

