

REF **AK-22123-F**

**ARROW**



**Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®, Arrowg+ard Blue® Catheter for High Volume Infusions**

**12** Fr. | **3** Lumen | **16** cm catheter length | **.035** inch dia. spring-wire guide

**Contents:**

- 1: Multi-Lumen Indwelling Catheter: 12 Fr. x 16 cm Radiopaque Polyurethane with Blue FlexTip®, Arrowg+ard® Antimicrobial Surface Treatment<sup>1</sup>, Extension Line Clamps
- 1: Spring-Wire Guide, Marked: .035" (0.89 mm) dia. x 23-5/8" (60 cm) (Straight Soft Tip on One End - "J" Tip on Other) with Arrow Advancer
- 1: Catheter: 18 Ga. x 2-1/2" (6.35 cm) Radiopaque over 20 Ga. RW Introducer Needle
- 1: Injection Needle: 25 Ga. x 1" (2.54 cm)
- 1: Blunt Fill Needle: 18 Ga. x 1-1/2" (3.81 cm)
- 1: Pressure Transduction Probe
- 1: Introducer Needle: 18 Ga. x 2-1/2" (6.35 cm) XTW and 5 mL Arrow® Raulerson Spring-Wire Introduction Syringe
- 1: Injection Needle: 22 Ga. x 1-1/2" (3.81 cm) and 5 mL Luer-Slip Syringe
- 1: Syringe: 3 mL Luer-Lock
- 1: Tissue Dilator: 12 Fr. (4.0 mm) x 14 cm
- 1: 5 mL 1% Lidocaine HCl Solution and Alcohol Prep
- 1: 3 mL Applicator 2% CHG and 70% IPA Chloraprep®<sup>2</sup> One-Step Solution with Hi-Lite Orange™ Tint (non-sterile solution)
- 1: SharpsAway® Disposal Cup
- 1: Drape: 24" x 36" (61 cm x 91 cm) with 4" (10 cm) fenestration
- 1: Scalpel: #11
- 2: Gauze Pad: 2" x 2" (5 cm x 5 cm)
- 5: Gauze Pad: 4" x 4" (10 cm x 10 cm)
- 1: Suture: 3-0 Braided Silk with Straight Needle


<sup>1</sup>Licensed under US Patent No. 6,706,024.  
<sup>2</sup>A registered trademark of CareFusion or one of its subsidiaries.

Rx only. **Warning: Read all package insert warnings, cautions, and instructions prior to use. Failure to do so may result in severe patient injury or death.**

**California Prop. 65**  
**Warning: Cancer and Reproductive Harm** [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Sterilized using ethylene oxide unless otherwise indicated in the contents list.  
 Fluid path components are non-pyrogenic.

**Contraindications: The Arrowg+ard Blue antimicrobial catheter is contraindicated for patients with known hypersensitivity to chlorhexidine, silver sulfadiazine and/or sulfa drugs (refer to product instructions for references).**



| Lumen             | Priming Volume* (mL) | Flow Rate† (mL/min) |
|-------------------|----------------------|---------------------|
| Distal (16 Ga.)   | 0.35                 | 51                  |
| Medial (12 Ga.)   | 1.20                 | 195                 |
| Proximal (12 Ga.) | 1.10                 | 254                 |

\* Priming volumes are approximate and are done without accessories.  
 † Based on 100 cm head height.



Not made with natural rubber latex



Store between 20-25°C (68-77°F).

Manufacturer:  
 Arrow International, Inc.  
 Subsidiary of Teleflex Incorporated  
 2400 Bernville Road  
 Reading, PA 19605 USA



www.teleflex.com  
 IK-22123-105A (2020-02)

REF **AK-22123-F**

**ARROW**



**Large-Bore Multi-Lumen Central Venous Catheterization Kit with Blue FlexTip®, Arrowg+ard Blue® Catheter for High Volume Infusions**

**12** Fr. | **3** Lumen | **16** cm catheter length | **.035** inch dia. spring-wire guide